

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Forcare B.V.	ForIndex	2.0	19-Sep-2008
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Audit Trail and Node Authentication (ATNA)	Secure Application	Encryption Option	
Consistent Time (CT)	Time Client		
Cross-Enterprise User Assertion (XUA)	X-Service-Provider		
Cross-Enterprise Document Sharing (XDS.a)	Document Registry		
Cross-Enterprise Document Sharing (XDS.b)	Document Registry	Patient Identity Feed	
Patient Demographics Query (PDQ)	Patient Demographics Supplier		
<b>Internet address for vendor's IHE information:</b> <a href="http://www.forcare.nl">http://www.forcare.nl</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	<a href="http://www.forcare.nl">http://www.forcare.nl</a>		
<b>DICOM</b>	<a href="http://www.forcare.nl">http://www.forcare.nl</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Forcare B.V.	ForStore	2.0	19-Sep-2008
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Audit Trail and Node Authentication (ATNA)	Secure Application	Encryption Option	
Consistent Time (CT)	Time Client		
Cross-Enterprise User Assertion (XUA)	X-Service-Provider		
Cross-Enterprise Document Sharing (XDS.a)	Document Repository		
Cross-Enterprise Document Sharing (XDS.b)	Document Repository		
<b>Internet address for vendor's IHE information:</b> <a href="http://www.forcare.nl">http://www.forcare.nl</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	<a href="http://www.forcare.nl">http://www.forcare.nl</a>		
<b>DICOM</b>	<a href="http://www.forcare.nl">http://www.forcare.nl</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Forcare B.V.	ForView	2.0	19-Sep-2008
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Audit Trail and Node Authentication (ATNA)	Secure Application	Encryption Option	
Consistent Time (CT)	Time Client		
Enterprise User Assertion (EUA)	Client Authentication Agent		
	Kerberized Server		
Cross-Enterprise User Authentication (XUA)	X-Service-User		
Cross-Enterprise Document Sharing (XDS.a)	Document Source		
	Document Consumer		
Cross-Enterprise Document Sharing (XDS.b)	Document Source		
	Document Consumer		
Cross-Enterprise Document Sharing for Imaging (XDS-I)	Imaging Document Consumer		
Patient Demographics Query (PDQ)	Patient Demographics Consumer		
Patient Identifier Cross-referencing (PIX)	Patient Identifier Cross-reference Consumer		
Scanned Documents (XDS-SD)	Content Consumer	View Option	
Sharing Laboratory Reports (XD-LAB)	Content Consumer	View Option	
Cross-Enterprise Sharing of Medical Summaries (XDS-MS)	Content Consumer	View Option	
<b>Internet address for vendor's IHE information: <a href="http://www.forcare.nl">http://www.forcare.nl</a></b>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	<a href="http://www.forcare.nl">http://www.forcare.nl</a>		
<b>DICOM</b>	<a href="http://www.forcare.nl">http://www.forcare.nl</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
Forcare B.V.	ForConnect	2.0	19-Sep-2008
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Audit Trail and Node Authentication (ATNA)	Secure Application	Encryption Option	
Consistent Time (CT)	Time Client		
Cross-Enterprise Document Sharing (XDS.a)	Document Source	Multiple Document Submission, Replace Document	
	Document Consumer		
Cross-Enterprise Document Sharing (XDS.b)	Document Source	Multiple Document Submission, Document Lifecycle Management	
	Document Consumer		
Cross-Enterprise Document Sharing for Imaging (XDS-I)	Imaging Document Source	Set of DICOM Instances PDF Report Text Report	
	Imaging Document Consumer		
Patient Identifier Cross-Referencing (PIX)	Patient Identifier Cross-Reference Consumer		
Scanned Documents (XDS-SD)	Content Creator		
Internet address for vendor's IHE information: <a href="http://www.forcare.nl">http://www.forcare.nl</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	<a href="http://www.forcare.nl">http://www.forcare.nl</a>		
<b>DICOM</b>	<a href="http://www.forcare.nl">http://www.forcare.nl</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	